



3008 Clarity

Specialty Table



Simplify Full Body Imaging and Patient Positioning



Make your imaging process effortless with the Skytron 3008 Clarity Specialty Table. The 3008 Clarity offers precise and infinite control of the tabletop's movement with its innovative floating top functionality along with 83 inches of radiolucent carbon fiber surface. With this wide range of motion and radiolucency, the Clarity is a versatile and adaptive table for a variety of specialties and case types.

Greater Imaging Access

- 83" fully radiolucent carbon fiber tabletop
- Allows clinicians to use the lowest doses for imaging and still maintain excellent image clarity
- 65" of unobstructed imaging window

Simple Dye Chasing

- Multi-directional slide on both X and Y axes
- The floating top function is engaged with a single button, allowing the user to reposition the top
- Reposition with ease even with patients up to 500lbs
- Ideal choice for vascular, orthopedic, and other imaging procedures




Multiple
Direction
Top-Slide



Unique
Procedures



Greater
Imaging
Access



Precise Hand Control

A specially designed, comfort grip hand control is all that is needed to articulate the 3008 Clarity floating top. The float function is engaged with a single button, allowing the user to push or pull the top in any direction with ease.

Precise Movement for Accurate Imaging

20° trend and tilt

Providing 20° of trend and tilt, the 3008 Clarity takes the work out of imaging setup. Tilt and trend support a 300 pound capacity through the entire range of motion.

Differentiating Features

- Radiolucent top for low dose images
- Fluoroscopy made easy with lateral and longitudinal movement of floating top
- One button float top control
- 83" Carbon-fiber top for full body imaging
- 65" of unobstructed imaging length
- One touch return to level top positioning
- 500 lb lift and top slide capacity



Skytron's Solutions



Sterile Processing

Processing Sinks and Stainless Offerings
Washers, Ultrasonics, and Cart Washers
Prep & Pack Tables & Sterilizers



Clinical

Surgical Tables
Booms & Lights
Full Offering of Stainless Products



Architectural & Environmental

Modular Walls
Configurable Window and Door Offerings
Built-in Cabinetry



Clinical Business Intelligence

Video Integration
Asset Tracking & Supply Chain Solutions
Surgical Displays

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV1 03-24
PN L2-010-156