



Sterling[®] C2000

Cart Washer-Disinfectors



Instrument Cart



Container Cart

Increase Productivity with Faster Process Times

The Sterling C2000 Cart Washer streamlines and automates the decontamination process – bringing unsurpassed efficiency and flexibility to the Sterile Processing Department. The C2000 is engineered to be environmentally friendly by conserving water, energy, and chemical consumption using its built-in control functions and water recirculation system.

Lower Cost of Ownership

- Lower water, chemical, and electric use than other washers
- Fast cycle times for increased throughput

Efficient Work Environment


- Noise level is less than 63 dB
- Optional steam-heated drying system

Maximize Your Space

- Four different sizes available lowers the need for custom construction
- Optional instrument cycles to augment existing washer-disinfectors



Water, Chemical,
& Electricity
Conservation



Low
Noise
Level



Fast
Process
Time

Efficient Water Coverage Throughout the Chamber

- 20 wash arms with 120 oscillating nozzles
- Pressure controlled based on load configuration
- Wide spray pattern for complete coverage
- Spray arms and nozzles can be cleaned without tools



The Sterling® Family of Washer-Disinfectors



The Sterling family of Washer-Disinfectors provides a full array of options for size and throughput to outfit any Sterile Processing Department with exactly what they need. Each model offers water, chemical, and electricity savings by design.

Comply With Your Medical Device IFU's

- Edit cycle parameters to meet your needs
- Document the performance of all cycles electronically

Boost Your Efficiency

- Best throughput per cycle in the least amount of space
- Choose the perfect sized solution for your location

Skytron's Solutions



Sterile Processing

Processing Sinks and Stainless Offerings
Washers, Ultrasonics, and Cart Washers
Prep & Pack Tables & Sterilizers



Clinical

Surgical Tables
Booms & Lights
Full Offering of Stainless Products



Architectural & Environmental

Modular Walls
Custom Window and Door Offering
UVC Disinfection Products



Clinical Business Intelligence

Video Integration
Asset Tracking & Supply Chain Solutions
Surgical Displays

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV2 09-22
PN L2-010-215