



Skytron

Disposable Camera Cover

Skytron's disposable camera covers are designed for use with Skytron surgical camera lights. The disposable covers are used over the camera and allow operating staff the ability to adjust the focus of the surgical light from within the sterile field. The covers come prepackaged and ready for use.

- Packaged 50 per case
- Container size - 21" W x 16" D x 19" H
- Product is sterilized using Ethylene Oxide
- Made in the USA
- Shelf life up to 5 years
- Latex free
- DEHP free
- BPA free
- No MSDS required



Skytron

Disposable Camera Cover



B1-715-65

- Simple snap-button engagement
- Single use for a Sterile Surface every procedure
- Skytron engineered to fit Stellar, Aurora LED, Aurora AUR, Aurora AUA, and Aurora AUT series light fixtures
- Clear lens area for uncompromised view and focus retention

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV2 09-22
PN L2-010-107