



Lumos

Surgical Light



Light is subjective. Lumos offers 60 different variations of light and 12 that produce 160,000 LUX. Lumos' personalized lighting enhances your ability to see what is important to you!

3 x 4 x 5

Spot
Sizes

Color
Temps

Adjustable
Intensities

Precisely Tailored Illumination™
Remembers You



Putting the Surgical Team in Control of Focus & Illumination

The Lumos surgical light is the latest LED innovation to the Skytron lighting portfolio. This next-generation surgical light puts enhanced clarity and focus in the clinician's hand. Utilizing our [best-in-class reflective LED technology](#), Lumos provides ideal illumination with [optimal shadow control](#).

Personalization

The Lumos surgical light brings crisp, deep-cavity illumination to any procedure with three spot sizes. With sixty combinations of spot size, color temperature, and intensity, this focusable light is sure to impress any surgeon. By utilizing the ten presets, Lumos creates efficiency and lowers set-up time and frustration.

Focusable

Controlling the focus from the sterile handle eliminates the need to reposition the lighthead up or down. The focal point adjusts according to the surgical team's preference — with only one twist of the sterile handle.

Intuitive Control

Lumos' flush-mounted wall control allows you to cycle through small, medium, and large spot sizes, 3800K, 4100K, 4500K, and 5000K color temperatures, five different intensities, and control the in-light surgical camera with a simple touch on the LCD screen.



Personalized
Lighting



Reflective LED
Technology



Intuitive
Control



Controlling spot size, color temperature, light intensity, and the camera from one intuitive control panel makes case setup easy and fast. Up to ten presets maximizes efficiency for the team.

Direct Focus: Less Movement, More Control



Greatest Intensity

Working Distance



Greatest Intensity

Working Distance

Adjust focus without moving the lighthead – one twist of the handle moves the focal point up or down and can be utilized at all 3 spot sizes.

Skytron's Solutions



Sterile Processing

Processing Sinks and Stainless Offerings
Washers, Ultrasonics, and Cart Washers
Prep & Pack Tables & Sterilizers



Clinical

Surgical Tables
Booms & Lights
Full Offering of Stainless Products



Architectural & Environmental

Modular Walls
Configurable Window and Door Offerings
Built-in Cabinetry



Clinical Business Intelligence

Video Integration
Asset Tracking & Supply Chain Solutions
Surgical Displays

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV2 03-24
PN L2-010-223