



Sterling P1000 & P1800

Washer-Disinfector Accessory Guide

Setting the standard for sterile processing

The Sterling P1000 and P1800 Washer-Disinfector series streamlines and automates the decontamination process – bringing unsurpassed efficiency and flexibility to the Sterile Processing Department.

Table of Contents	1
Wash Racks.....	1
Rack Inserts and Support.....	5
Injectors and Nozzles.....	7
P1000 Specialty Accessories.....	12
P1800 Specialty Accessories.....	13

Wash Racks

Flexible Rack

Three Levels

- Wash in up to three levels
- Includes grids and detachable washing arms
- Maximum height of 5.1 inches (13 cm) on the top level and 7.8 inches (20 cm) on each lower level

Item #1000-FLX3

Item #1800-FLX3



Flexible Rack

Four Levels

- Wash in up to four levels
- Includes grids and detachable washing arms
- Maximum height of 4.6 inches (11.8 cm) on the top level and 4.6 inches (11.8 cm) on each lower level

Item #1000-FLX4

Item #1800-FLX4



Flexible Rack

Five Levels

- Wash in up to five levels
- Includes grids and detachable washing arms
- Can wash up to 15 DIN instrument baskets per load
- Maximum height of 3.2 inches (8.2 cm) between levels

Item #1000-FLX5

Item #1800-FLX5



Wash Racks

Fixed Rack

Four Levels

- Washes up to 12 DIN instrument baskets per load
- The loading space in each level is up to 3.2 inches (8.2 cm)

Item #1000-FIX4

Item #1800-FIX4



Fixed Rack

Five Levels

- Washes up to 15 DIN instrument baskets per load
- The loading space in each level is up to 2.5 inches (6.3 cm)

Item #1000-FIX5

Item #1800-FIX5



Fixed Rack

Six Levels

- Washes up to 18 DIN instrument baskets per load
- Maximum height of 2 inches (5 cm) at top, max 2.5 inches (6.3cm) all lower

Item #1000-FIX6

Item #1800-FIX6



Container Rack

The vertical container load carrier is designed to wash a number of containers within the wash rack footprint.

For medium containers sized 600 x 300 x 300 mm

Item #1000-CNTR-MED

Item #1800-CNTR-MED



Wash Racks

Container Rack

The vertical container load carrier is designed to wash a number of containers within the wash rack footprint.

For small containers sized 600 x 300 x 150 mm

Item #1000-CNTR-SM

Item #1800-CNTR-SM



Insert Rack

Baby bottle rack

Insert for efficiently washing up to 96 baby bottles per load.

- 125 ml bottles
- Flexible and removable insert
- Used in conjunction with two level wash racks

Item #1000-BB-125




Insert Rack

Baby bottle rack

Insert for efficiently washing up to 96 baby bottles per load.

- 250 ml bottles
- Flexible and removable insert
- Used in conjunction with two level wash racks

Item #1000-BB-250



Insert Rack

Washing OR Shoes

Insert for efficiently washing up to 16 OR shoes per load.

- 16 pairs of hospital clogs or shoes
- Flexible and removable insert

Item #1000-FLX-SH

Item #1800-FLX-SH



Wash Racks

MIC Rack

58 Basic Nozzles

MIC rack complete, grommet, with 56 basic nozzles.

Item #1000-MIC-WN

Item #1800-MIC-WN



MIC Rack

Without Nozzles

MIC rack without nozzles, and grommet.

Item #1000-MIC-WON

Item #1800-MIC-WON



Da Vinci Rack

Da Vinci rack without holders Si, Xi.

The example rack is fitted with various holders and a RIFD Reader. The capacity of this rack is 14 holders.

Four containers

Item #1000-DVXISI

Item #1800-DVXISI



Anesthesia Rack

Tubes and Equipment

Designed to wash anesthesia tubes and equipment. Rack holds up to 12 anesthesia tubes, breathing bags, face masks, intubation instruments, laryngoscopes etc. One level suitable for DIN basket for face masks, intubation instruments, laryngoscopes and other equipment.

Item #1000-AN

Item #1800-AN



Rack Inserts and Support

Shelf Level Rack

Washing Stainless Bowls

Shelf level for Sterling P1000 and P1800 wash racks for washing and disinfecting kidney bowls.

- Capable of holding up to six bowls
- Can hold bowls up to a max diameter of 13.8 inches (35 cm)

Item #1000-FLX-KB



Shelf Level Rack

Washing Stainless Bowls

Support for four washing bowls three per level for five and four level flexible racks.

- Four washing bowls
- Washing capacity of 12 bowls of different sizes
- Maximum diameter of 14.2 inches (36 cm) per bowl

Item #1000-FLX-BWL



Spike For Respirator

Anesthesia Rack

Spike for respirator bellow for Anesthesia rack.

Item #2908265



Da Vinci Si Holder

Da Vinci Si holder 8/5 mm 10 per rack

14 holders per rack

Item #1000-SI85



Rack Inserts and Support

Da Vinci Xi

Da Vinci Xi holder 8 mm and single site instrument.

14 holder per rack

Item #1000-XI8



Da Vinci Xi

Stapler holder for Da Vinci Xi.

Ten holders per rack.

Item #1000-XIST



Da Vinci Single Site

14 holders per rack

Item #1000-SS



Da Vinci Si

Da Vinci Si stapler holder

14 holders per rack

Item #1000-SIST




Rack Inserts and Support

Fine Mesh Grate For Racks

P1000 Item #2908469

P1800 Item #2918469



Cleaning Kit

Kit for cleaning power tools.

Item #2908314



Injectors and Nozzles

Fine Filter Module

For Dental and Ophthal Units

Fine filter module for Dental and Ophthal units without units for 5 level racks.

P1000 Item #2908539


P1800 Item # 2918539



Dental Unit 016/020

- Four - 016
- Four - 020
- Three units per module

Item #2908508



Injectors and Nozzles

Dental Unit 016

- Eight piece 016 adapter
- Three units per module

Item #2908432




Distributor Unit

Adapters

- Eight piece 020 adapters
- Two unit per module.

Item #2908506



Injector Spikes 04

Anesthesia Rack

Injector spikes for tubes 04 for Anesthesia rack.

Item #2908255



Dental Unit 020

- Four - 016
- Four - 020
- Three - units per module

Item #2908426



Injectors and Nozzles

Nozzle 04

MIC Rack

Injector nozzle screw 04 for old MIC rack and module with screw.

Item #2908364



Injector Nozzle 02.6

MIC Rack

Injector nozzle 02.6 (grommet) length 50 mm. for MIC rack and modules.

Item #2908441



Injector Nozzle Closed

MIC Rack

Injector nozzle 04 closed top (grommet) length 120mm. for MIC rack and modules.

Item #2908426



Injector Nozzle Open

MIC Rack

Injector nozzle 04 open top (grommet) length 120mm. for MIC rack and modules.

Item #2908428



Injectors and Nozzles

Basic Nozzle

MIC Rack

Basic nozzle screw for tubular instruments for old MIC rack and module with screw.

Dimensions 10 x 6 x 2 mm.

Item #2900469



Basic Nozzle - Long

MIC Rack

Basic nozzle, long, (grommet) for tubular instruments.

- Maximum O13 mm.
- Sleeve length 55 mm.

Item #2900807



Basic Nozzle - Short

MIC Rack

Basic nozzle, short, (grommet) for tubular instruments.

- Maximum O13 mm.
- Sleeve length 27 mm.

Item #2900808



Ophthal Nozzle

Ophthal nozzle fine mesh filter

Item #2908373



Injectors and Nozzles

Caps For Basic Nozzles


1: Caps for basic nozzle 02 - Item #2900804

2: Caps for basic nozzle 03 - Item #2900805

3: Caps for basic nozzle 04 - Item #2900806

4: Caps for basic nozzle 03-08 - Item #2900802

5: Caps for basic nozzle 06-011 - Item #2900803



Universal Nozzle

MIC Rack

Nozzle universal screw for tubular instruments for old MIC rack and module with screw.

Item #2908184



Dental Nozzle

Dental nozzle 016


Item #2908303



Dental Nozzle 020

Dental nozzle 020

Item #2908304



Nozzle 04 Closed

MIC Rack

Injector nozzle screw 04 closed for old MIC rack and module with screw.

Item #2908283



Injectors and Nozzles

Nozzle 02.6

MIC Rack

Injector nozzle screw 02.6 for old MIC rack and module with screw.

Item #2908439



Ophthal Unit

Three units per module

Item #2908519



P1000 Specialty Accessories

Manual Transport Cart

Manual loading cart transports Sterling P1000 wash racks in all configurations. Easily maneuverable around all environments.

- Fits all Skytron Sterling P1000 wash rack configurations
- Built in drip tray and collector for excess water
- Easy to use coupling system for connection to the Sterling P1000 washer

Item #1000-MC



P1800 Specialty Accessories

Flexible Rack

5 Levels

Universal three level rack suitable for various OR equipment from orthopedic surgeries - like hip, foot, knee, hand, shoulder or other surgeries that involve numerous instrument sets.

Item #1800-FLX5-PLUS



Da Vinci Xi Optics

Holder / Scope Holder

Da Vinci Xi Optics Holder/ Scope Holder.

Four holders

Item #1800-DVXI



Manual Transport Cart

Manual loading cart transports Sterling P1800 wash racks in all configurations. Easily maneuverable around all environments.

- Fits all Skytron Sterling P1800 wash rack configurations
- Built in drip tray and collector for excess water
- Easy to use coupling system for connection to the Sterling P1800 washer

Item #1800-MC



P1800 Specialty Accessories

Automatic Loading Cart

The Automatic loading system for the Sterling P1800 can be connected to the clean and unclean side of the machine.

- Fits all Skytron Sterling P1800 wash rack configurations
- Built in drip tray and collector for excess water
- Easy to use coupling system for connection to the Sterling P1800 washer
- Minimize heavy handling of instruments

Item #1800-ALT



Automatic Unloading Cart

The Automatic unloading system for the Sterling P1800 can be connected to the clean and unclean side of the machine.

- Fits all Skytron Sterling P1800 wash rack configurations
- Built in drip tray and collector for excess water
- Easy to use coupling system for connection to the Sterling P1800 washer
- Minimize heavy handling of instruments

Item #1800-AULT



Skytron's Solutions



Sterile Processing

Processing Sinks and Stainless Offerings
Washers, Ultrasonics, and Cart Washers
Prep & Pack Tables & Sterilizers



Clinical

Surgical Tables
Booms & Lights
Full Offering of Stainless Products



Architectural & Environmental

Modular Walls
Custom Window and Door Offering
UVC Disinfection Products



Clinical Business Intelligence

Video Integration
Asset Tracking & Supply Chain Solutions
Surgical Displays

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV2 09-22
PN L2-010-198