

# Sterling P1800

Washer-Disinfector





## Process your medical devices the most efficient, and compliant way possible

The Sterling P1800 is designed to set a new standard for Sterile Processing Departments. With a large capacity, small footprint, fast cycles, low water and electrical consumption and most important, disinfected devices on a consistent basis.

### Lower Cost of Ownership

- Less water and chemical use than other washers
- Save 50% in chemical and sewage costs

### Comply With Your Medical Device IFU's

- Edit cycle parameters to meet your needs
- Document the performance of all cycles electronically

### Maximize Your Efficiency


- Best throughput per cycle in the least amount of space
- Up to 18 DIN trays per cycle



Water  
Conservation



Customizable  
Cycles



Conserve  
Space



**Automate your SPD** with RFID equipped manifold racks that select and start your next cycle.

**Accommodate instrument trays of any size** P1800's modular manifold racks are easily configured.

**Maintain consistent performance** with a daily on-screen user checklist and remote service connectivity.

**Protects against water quality changes** utilizing an optional conductivity sensor.

**Elevate your disinfection process** with optional preheat boosters to shorten cycles and recirculation tanks to save critical water.

## Configured for Versatility

Instruments



Utensils



MIS Specific



# Skytron's Solutions



## Sterile Processing

Processing Sinks and Stainless Offerings  
Washers, Ultrasonics, and Cart Washers  
Prep & Pack Tables & Sterilizers



## Clinical

Surgical Tables  
Booms & Lights  
Full Offering of Stainless Products



## Architectural & Environmental

Modular Walls  
Custom Window and Door Offering  
UVC Disinfection Products



## Clinical Business Intelligence

Video Integration  
Asset Tracking & Supply Chain Solutions  
Surgical Displays

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

**To learn more, visit our website at [www.skytron.com](http://www.skytron.com)**



5085 Corporate Exchange Blvd. SE  
Grand Rapids, MI 49512  
1.800.SKYTRON (759-8766)  
email: [info@skytron.com](mailto:info@skytron.com)  
[www.skytron.com](http://www.skytron.com)



REV1 09-22  
PN L2-010-196