



Processing Sink

Height-Adjustable Solution for Sterile
Processing Departments

Ergonomically friendly. Sanitary. Efficient.

Skytron's Processing Sinks complement your disinfection processes and provide a means for the pre-cleaning, inspection and sorting of surgical instruments.

Standard Features

- Easy access to connection points with open base
- Single, dual, or triple basins
- Fully welded construction with reinforced tops
- Right to left workflow or left to right workflow
- 37" standard height for stationary sinks
- 4" backsplash
- Under shelf
- Adjustable feet
- Etched gallon markings
- Basin size: 24"W x 10"D x 17"L

Optional Features

- Ergonomically friendly electronic adjustable height range: 35"–42"
- Mixing valve faucet
- Pre-rinse faucet
- DI water faucet
- Air gun kit with shut off valve and 5 tips
- Water gun kit with shut off valve and 8 tips
- Drain boards
- 6", 8", or 12" backsplash
- 4", 6", 8", or 12" side splash
- Hand lever drains
- Foot operated drains (not available for adjustable height units)
- Copper drain system (not compatible with deionized water)
- 30"W x 10"D x 17"L sink upgrade



Alleviate Physical Stress

Height adjustable option ranging from 35"-42"



Cleaner Work Environment

With laser etched gallon markings, there are fewer areas for bacteria to hide

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV3 01-24
PN L2-010-178