



Stainless

Scrub Sinks

Skytron's Stainless Steel Scrub Sinks provide quality, convenience and flexibility for every clinical need. Scrub Sinks are available with knee-kick and/or optional infrared modules for "hands-free" operation. In addition, the optional Eye Wash (EW) feature provides a central, visible location to all staff and saves valuable floor and wall space.



Triple Bay

MK Sink: 96"W 27-1/2"D 37"H

ES Sink: 76"W 20-1/2"D 23-9/16"H



Dual Bay

MK Sink: 64"W 27-1/2"D 37"H

ES Sink: 47"W 20-1/2"D 23-9/16"H



Single Bay

MK Sink: 32"W 27-1/2"D 37"H

ES Sink: 25"W 20-1/2"D 23-9/16"H

Quick Ship Option

SS2121-MK-SKY: 32"W 27-1/2"D 37"H

Ready-to-ship sink with infrared control, eyewash, timer, dual laminar flow, blending valve, thermostatic mixing valve, and a dual bay overhead wall shelf

Quick Ship Option

SS2121-MK-SKY: 32"W 27-1/2"D 37"H

Ready to ship sink with infrared control, eyewash, timer, dual laminar flow, blending valve, and a single bay overhead wall shelf



Stainless Scrub Sinks

Features	Standard		Optional
	MK	ES	MK
Single, dual, or triple stations	X	X	
Deep sloping basin	X	X	
Adjustable thermostatic controlled pressure regulating mixing valve (one per station)	X		
Has a preset mixing valve		X	
Knee operated water	X	X	
Check mixing valve stop strainers for ease of maintenance	X		
High mount swivel gooseneck spout with adjustable swivel spray aerator	X	X	
Flow control valve	X	X	
Backflow preventers	X	X	
Infrared sinks are equipped with a 24 volt plug-in transformer	X	X	
Splash screen dividers with dual and triple station sinks	X		
Main support and Z-brackets	X	X	
304 stainless steel polished and ground to a #4 finish	X	X	
Infrared operation			X
Knee operated water with IR water operation			X
Knee operated soap with IR water operation			X
Foot operated water			X
Foot operated soap			X
Digital scrub timer			X
Eye wash			X
Thermostatic mixing valve			X
Add-on to Eye Wash option			X
Laminar flow			X
220/240V power supply			X
Pedestal base			X
In-wall chair carrier			X
Overhead wall shelf over sink			X

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
 Grand Rapids, MI 49512
 1.800.SKYTRON (759-8766)
 email: info@skytron.com
www.skytron.com



REV3 10-23
 PN L2-010-218