



Stainless

Prep & Pack Tables



Ergonomic Design for Efficient Assembly & Packaging in Your SPD

Skytron Prep & Pack tables are designed to alleviate the physical stress placed on technicians performing the inspection, assembly, and packaging in your Sterile Processing Department. Standing in one position for hours on end, or using a workstation designed for someone else, can lead to uncomfortable neck or back pain and unsatisfied employees. Put your faculty's well-being in mind, and choose a Prep & Pack table that easily adjusts to the various physical attributes and requirements of the human body.

Specifications	Prep & Pack TWS-3060	Prep & Pack TWS-3072	Prep & Pack TWS-3660	Prep & Pack TWS-3672
Frame Height	84" - 92"	84" - 92"	84" - 92"	84" - 92"
Width	60"	72"	60"	72"
Depth	32.25"	32.25"	38.25"	38.25"
Desk Height	36" - 44"	36" - 44"	36" - 44"	36" - 44"
Height Adjustment Options	Fixed, Manual, and Powered	Fixed, Manual, and Powered	Fixed, Manual, and Powered	Fixed, Manual, and Powered



Comfortable Workstation



Durable Stainless Steel



Easy Accessory Attachment

Differentiating Features

- 36"-44" height adjustable work surface
- 43" vertical riser frame raises and lowers with the base
- Fixed height, manual, powered models available
- Ergonomically correct alignment of computer monitors, lighting, bins, and pegboards
- Heavy duty inner/outer leg extrusions for increased stability
- 500lb dynamic and static weight capacity



Optional Accessories



H0501 Monitor Mount




H0509-03 Magnifier



TWS-ASH-B Angled Shelf



TWS-BR60-B Bin Rail



TWS-CPU-30 CPU Holder



TWS-DWR Drawer



TWS-OSB60-B Power Strip



TWS-KT Keyboard Tray



TWS-OLT-60 Light Bar



TWS-LS3060-B Lower Shelf



TWS-PB-60 Pegboard



TWS-RH60-B Spool



TWS-SH-60-B Shelf Adjustment



TWS-WB-60 Whiteboard



IC0128-01
3-3/32" SS single peg board hook

IC0128-02

5-3/16" SS single peg board hook

IC0128-03

8-3/4" SS single peg board hook

IC0128-05

9" SS single peg board hook

IC0128-09

9" SS dual peg board hook

IC0901

8" SS single peg board hook

Skytron's Solutions



Sterile Processing

Processing Sinks and Stainless Offerings
Washers, Ultrasonics, and Cart Washers
Prep & Pack Tables & Sterilizers



Clinical

Surgical Tables
Booms & Lights
Full Offering of Stainless Products



Architectural & Environmental

Modular Walls
Configurable Window and Door Offerings
Built-in Cabinetry



Clinical Business Intelligence

Video Integration
Asset Tracking & Supply Chain Solutions
Surgical Displays

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service.

Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.

To learn more, visit our website at www.skytron.com



5085 Corporate Exchange Blvd. SE
Grand Rapids, MI 49512
1.800.SKYTRON (759-8766)
email: info@skytron.com
www.skytron.com



REV2 03-24
PN L2-010-221