



# Skytron

Warming Cabinets

# Long Life. High Quality. Accurate Temperatures.



Accelerate patient recovery time and reduce surgical site infections by maintaining normothermia from admission to discharge. Designed to thoroughly heat and maintain IV fluids, irrigation fluids and blankets to healthcare or manufacturer temperature guidelines, Skytron's broad range of warming cabinet configurations also meet storage, space, and application requirements throughout a facility. Simplify tasks and enhance convenience with choices of remote temperature monitoring, touch screen control or glass doors for visualization.

## Safe and Secure

- LED temperature display with keyed lockout and data logging
- Monitor, control, report, and troubleshoot all from a centralized location (touch screen only)
- Touch screen option allows the user to easily control the temperature permissions

## Better Patient Outcomes

- Maintain accurate temperatures with our robust magnetic closure system
- Perforated shelves enhance airflow to ensure contents are equally warmed

## Low Cost of Ownership

- Made from durable 300 series stainless steel throughout
- Premium components assure optimal performance for the life cycle of the warmer



Safe and  
Secure



Accurate  
Temperatures



Low  
Cost of  
Ownership

# Single Chamber Units



Product Number	Overall Dimensions	Approximate Capacity
W1-252421	24.5"H x 24" W x 20.5" D	9 Blankets or 8 Solution Bottles (1 Liter Each)
W1-252427	24.5"H x 24" W x 26.5" D	9 Blankets or 12 Solution Bottles (1 Liter Each)
W1-253021	24.5"H x 30" W x 20.5" D	18 Blankets or 8 Solution Bottles (1 Liter Each)
W1-253027	24.5"H x 30" W x 26.5" D	18 Blankets or 20 Solution Bottles (1 Liter Each)
W1-362421	36"H x 24" W x 20.5" D	16 Blankets or 20 Solution Bottles (1 Liter Each)
W1-362427	36"H x 24" W x 26.5" D	16 Blankets or 25 Solution Bottles (1 Liter Each)
W1-363021	36"H x 30" W x 20.5" D	32 Blankets or 36 Solution Bottles (1 Liter Each)
W1-363027	36"H x 30" W x 26.5" D	32 Blankets or 40 Solution Bottles (1 Liter Each)
W1-652421	64.75"H x 24" W x 20.5" D	32 Blankets or 48 Solution Bottles (1 Liter Each)
W1-652427	64.75"H x 24" W x 26.5" D	80 Blankets or 120 Solution Bottles (1 Liter Each)
W1-653021	64.75"H x 30" W x 20.5" D	66 Blankets or 72 Solution Bottles (1 Liter Each)
W1-653027	64.75"H x 30" W x 26.5" D	80 Blankets or 120 Solution Bottles (1 Liter Each)
W1-752421	74.5"H x 24" W x 20.5" D	32 Blankets or 60 Solution Bottles (1 Liter Each)
W1-752427	74.5"H x 24" W x 26.5" D	80 Blankets or 120 Solution Bottles (1 Liter Each)
W1-753021	74.5"H x 30" W x 20.5" D	80 Blankets or 72 Solution Bottles (1 Liter Each)
W1-753027	74.5"H x 30" W x 26.5" D	80 Blankets or 120 Solution Bottles (1 Liter Each)

# Dual Chamber Units



Product Number	Overall Dimensions	Approximate Capacity
W2-652421	64.75"H x 24" W x 20.5" D	26 Blankets or 36 Solution Bottles (1 Liter Each)
W2-652427	64.75"H x 24" W x 26.5" D	52 Blankets or 80 Solution Bottles (1 Liter Each)
W2-653021	64.75"H x 30" W x 20.5" D	50 Blankets or 54 Solution Bottles (1 Liter Each)
W2-653027	64.75"H x 30" W x 26.5" D	52 Blankets or 96 Solution Bottles (1 Liter Each)
W2-752421	74.5"H x 24" W x 20.5" D	32 Blankets or 48 Solution Bottles (1 Liter Each)
W2-752427	74.5"H x 24" W x 26.5" D	52 Blankets or 80 Solution Bottles (1 Liter Each)
W2-753021	74.5"H x 30" W x 20.5" D	52 Blankets or 72 Solution Bottles (1 Liter Each)
W2-753027	74.5"H x 30" W x 26.5" D	52 Blankets or 120 Solution Bottles (1 Liter Each)

# Skytron's Solutions



## Sterile Processing

Processing Sinks and Stainless Offerings  
Washers, Ultrasonics, and Cart Washers  
Prep & Pack Tables & Sterilizers



## Clinical

Surgical Tables  
Booms & Lights  
Full Offering of Stainless Products



## Architectural & Environmental

Modular Walls  
Configurable Window and Door Offerings  
Built-in Cabinetry



## Clinical Business Intelligence

Video Integration  
Asset Tracking & Supply Chain Solutions  
Surgical Displays

Skytron provides safe and efficient healthcare solutions delivered through a local representative network that values long-term relationships and comprehensive service. Focused on providing the highest quality solutions while delivering the lowest cost of ownership, we prioritize transparency and accountability in our actions.  
**To learn more, visit our website at [www.skytron.com](http://www.skytron.com)**



5085 Corporate Exchange Blvd. SE  
Grand Rapids, MI 49512  
1.800.SKYTRON (759-8766)  
email: [info@skytron.com](mailto:info@skytron.com)  
[www.skytron.com](http://www.skytron.com)

